






Vacuum Manifolds for 96-Well Plates

The Orochem Vacuum Manifold is an anodized aluminum manifold that has been designed for vacuum filtration of the OmniSorb filter plates, Orochem 2000 filter plates and removable well format plates.

- Replaceable Gasket –Ensures perfect vacuum seal. The gasket can be replaced instead of entire manifold.
- Vacuum Gauge provides precise reading of vacuum level.
- Chemically Resistant Anodized Base Prolongs life of product.
- Compatibility with most manufacturers' Plates – SBS Compatible footprint.
- Anti Slip Base Prevents slippage and provides improved safety during operation.
- Small footprint ensures minimal space requirement on automation instrument deck
- Additional Vacuum Control options available



1 ml 96 well filter plate



2 ml 96 well filter plate

Orochem Technologies Inc.
331, S. Eisenhower Ln, Lombard, IL 60148
Phone: 630-916-0225 Fax: 630-916-0250
Email: sales@orochem.com